


Gluteus Implants

by EUROSILICONE®

ES 62
Smooth Surface, Round




Product Code	Diameter (cm)	Projection (cm)	Volume (cc)
62140	8.4	4.1	140
62180	9.1	4.3	180
62200	9.4	4.4	200
62220	9.7	4.5	220
62240	10.3	4.7	240
62260	10.4	4.7	260
62280	10.6	4.9	280
62300	10.8	5.0	300



Ürüne Genel Bakış

- 2 farklı konfigürasyon.
- Yuvarlak ve anatomik form.
- Yüksek kohesivite.
- Pürüzsüz yüzey.

ES 60
Smooth Surface, Anatomical



Product Code	Height (cm)	Width (cm)	Projection (cm)	Volume (cc)
60150	13.1	9.6	2.3	150
60210	14.6	10.7	2.6	210
60230	12.0	11.2	3.5	230
60275	12.4	11.4	3.8	275
60300	13.9	10.9	3.6	300
60325	14.7	11.3	3.7	325
60375	14.7	11.7	4.1	375
60425	15.0	12.4	4.5	425
60480	15.5	13.2	4.8	480

