


# Remote Valve Tissue Expanders

by EUROSILICONE®

## Ürüne Genel Bakış

- Yuvarlak form.
- Dikdörtgen form.
- Silindirik kavisli.
- Silindirik düz.
- Hilal form.



ES 12  
Round Form

| Product Code | Diameter (cm) | Height (cm) | Valve (cm) | Volume (cc) |
|--------------|---------------|-------------|------------|-------------|
| 12010        | 3.1           | 2.1         | 2.0        | 10          |
| 12120        | 7.1           | 4.8         | 2.0        | 120         |
| 12150        | 9.5           | 4.0         | 2.0        | 150         |
| 12200        | 7.9           | 5.1         | 2.8        | 200         |
| 12250        | 9.0           | 5.7         | 2.8        | 250         |
| 12300        | 10.0          | 6.3         | 2.8        | 300         |
| 12375        | 13.0          | 5.0         | 2.8        | 375         |
| 12400        | 11.0          | 6.0         | 2.8        | 400         |
| 12500        | 12.8          | 6.3         | 2.8        | 500         |
| 12600        | 13.8          | 6.7         | 2.8        | 600         |
| 12700        | 14.8          | 7.5         | 2.8        | 700         |
| 12800        | 15.7          | 6.7         | 2.8        | 800         |
| 12900        | 16.6          | 6.5         | 2.8        | 900         |
| 121100       | 17.4          | 8.2         | 2.8        | 1100        |


All device dimensions correspond to the suggested fill volume.

GCA™ Remote Valve Tissue Expanders must be filled with a sterile isotonic solution only.

A fill tube is connected to every tissue expander to allow inflation of the device through an injection port.


All GCA™ Remote Valve Tissue Expanders with a recommended fill volume of more than 150cc are provided with a 28mm injection port, smaller volume devices are provided with a 20mm injection port.

ES 14  
Rectangular Form




| Product Code | Length (cm) | Width (cm) | Projection (cm) | Valve (cm) | Volume (cc) |
|--------------|-------------|------------|-----------------|------------|-------------|
| 14016        | 4.0         | 2.1        | 2.1             | 2.0        | 16          |
| 14035        | 4.0         | 2.6        | 3.6             | 2.0        | 35          |
| 14050        | 5.0         | 3.0        | 3.6             | 2.0        | 50          |
| 14100        | 4.9         | 3.9        | 5.1             | 2.0        | 100         |
| 14150        | 5.8         | 4.8        | 5.0             | 2.0        | 150         |
| 14200        | 8.5         | 5.0        | 4.9             | 2.8        | 200         |
| 14300        | 10.0        | 6.0        | 6.3             | 2.8        | 300         |
| 14400        | 10.8        | 5.8        | 5.9             | 2.8        | 400         |
| 14520        | 11.7        | 7.8        | 5.4             | 2.8        | 520         |
| 14650        | 11.7        | 9.0        | 5.8             | 2.8        | 650         |
| 14720        | 14.5        | 7.8        | 5.8             | 2.8        | 720         |
| 14750        | 13.5        | 8.7        | 5.9             | 2.8        | 750         |
| 14800        | 15.5        | 8.2        | 5.8             | 2.8        | 800         |
| 14950        | 17.5        | 8.8        | 5.5             | 2.8        | 950         |

ES 16  
Cylindrical Straight




| Product Code | Length (cm) | Width (cm) | Projection (cm) | Valve (cm) | Volume (cc) |
|--------------|-------------|------------|-----------------|------------|-------------|
| 16060        | 10.0        | 3.9        | 2.7             | 2.0        | 60          |
| 16100        | 11.0        | 3.5        | 3.5             | 2.0        | 100         |
| 16150        | 10.0        | 5.4        | 4.0             | 2.0        | 150         |
| 16200        | 14.0        | 4.5        | 4.5             | 2.8        | 200         |
| 16250        | 10.2        | 7.2        | 5.4             | 2.8        | 250         |
| 16300        | 16.0        | 5.0        | 5.0             | 2.8        | 300         |
| 16400        | 17.5        | 6.0        | 6.0             | 2.8        | 400         |

ES 16C  
Cylindrical Curved



| Product Code | Length (cm) | Width (cm) | Projection (cm) | Valve (cm) | Volume (cc) |
|--------------|-------------|------------|-----------------|------------|-------------|
| 16C250       | 14.5        | 7.7        | 4.5             | 2.8        | 250         |

ES 19  
Croissant Form



| Product Code | Length (cm) | Width (cm) | Projection (cm) | Valve (cm) | Volume (cc) |
|--------------|-------------|------------|-----------------|------------|-------------|
| 19060        | 8.5         | 6.0        | 3.9             | 2.0        | 60          |
| 19100        | 10.8        | 8.5        | 5.5             | 2.0        | 100         |
| 19125        | 10.0        | 7.0        | 4.8             | 2.0        | 125         |
| 19250        | 13.0        | 10.0       | 5.7             | 2.8        | 250         |
| 19500        | 17.0        | 12.0       | 7.8             | 2.8        | 500         |